## CB5TR/L

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*The crown-cutters CB5TR and CB5TRL are distinguished by their impressive versatility – they are capable of efficiently separating both metal crowns and crowns with ceramic veneers in record time. The instruments are at their most effective when used for cutting a crown in several steps. Apply the crown cutter several times to cut off smaller sections with amazing efficiency.*

*These crown cutters can be identified by the blue ring on their shank. The instruments are available in the sizes 012 and 014.*

*Our recommendation:  
Provided with a longer working part of 5 mm, the CB5TRL is the latest addition to our range. This new version is ideally suitable for separating larger crowns in the lateral region.*

### Advantages:

- Special diamond toothing
- The tapered working part adapts particularly well to the surface
- Crowns, bridges and veneers can be cut quickly and with ease
- Suitable for all commonly used metal alloys and veneers made of low-fusion ceramics and titanium

### Recommendations for use:

- Suitable for micro-motors and dental turbines
- Optimum speed in both power systems:  $\varnothing_{opt.} 160,000 \text{ rpm}$
- Always work with cooling spray (at least 50 ml/min.)
- Moderate contact pressure of 2N
- For cutting ceramic veneers or metal crowns, the instrument should be applied to the surface at an angle of approx. 45°, i.e. the crown-cutter hits the surface to be cut in a tilted position.

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
*The All-Rounder  
Multifunktional  
for ceramic veneers  
and metal alloys*

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
*The Turbo  
for non-precious metal  
alloys and gold*

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
*Aggressive –  
for non-precious metal  
alloys and gold*

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
*The Multifunctional Eco  
for non-precious metal  
alloys and gold*

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DIAMANT GmbH