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### Two-step spiral polishers for polishing composites

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EN

The two new spiral polishing wheels P9480M/F provided by D+Z are intended for polishing composite fillings after shaping with finishing instruments. They create a natural looking high shine finish in just two polishing steps.

- Step 1: Pre-polish with the dark grey polishing wheel P9480M
- Step 2: High-shine finish with the light grey polishing wheel P9480F

Thanks to their special shape, these spiral polishing wheels are highly flexible and can, therefore, be used universally. The polishers adapt to different surface structures so that all areas are reached with ease: occlusal, incisal, labial, buccal, lingual, proximal.


These polishing wheels are intended for single use and are unsuitable for reprocessing.

#### Advantages:

- Flexible spiral shape with 22 lamellae
- High shine in just two steps
- Integrated diamond abrasive particles
- Suitable for all current composites

#### Recommendations for use:

- Pre-polishing: dark grey
- High-shine polishing: light grey
- To be used in the micromotor without polishing paste
- Optimum speed:  $\varnothing_{opt}$  6,000 rpm
- Apply cooling agent
- Apply minimum contract pressure. The lamellae should not be pressed together completely
- Single-use product

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© D+Z · Germany · 03/19/419507V1
DIAMANT GmbH