Ausborenen alter Füllungen, Removal of old fillings,
Retrait des obturations, Remoción de obturaciones,
Escavazione di vecchie otturazioni

Verpackungseinheit, Packing unit, Unité d'emballage,
Unidad de embalaje, Unità d'imballo

Enthält gefährliche Substanzen, Contains hazardous
substances, Contient des substances dangereuses,
Contiene sustancias peligrosas, Contiene sostanze
pericolose

CAS 7440-48-4: Cobalt, 7440-48-4: Cobalt, 7440-48-4:
Cobalt, 7440-48-4: Cobalt, 7440-48-4: Cobalt


Fig. 1


Fig. 2

Amalgamentferner

1. Indikation

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Tungsten carbide instrument for removing amalgam fillings.

2. Contraindication

No product-related contraindications are known.

3. Start-up

- Those medical devices that are supplied non-sterile have to be prepared prior to first use.
- The instrument can be used in the turbine or micromotor.
- Heat generation and development of mercury vapor is largely reduced when used in the micromotor at a speed of 160,000 rpm.
- Apart from this general information, please observe the individual speed indicated on the packaging.
- The maximum permissible speed indicated on the package must not be exceeded in order to avoid risk of injury due to instrument breakage.
- Use low contact pressure.
- Never allow the instruments to jam. Do not use instruments for leverage.
- Drill with care and control to avoid loss of sound tooth substance as the instrument has a highly efficient blade configuration.
- Make sure that there is always sufficient coolant available. A minimum amount of cooling water of 50 ml/min is necessary to avoid damage to the pulp and increased mercury vapour.

4. Proper use

- Drilling is carried out step by step:
- Drill into the amalgam filling in an axial direction or in an inclined position (Fig. 1).
- Then cut several longitudinal and transversal separation grooves, depending on the size of the filling (Fig. 2). By dividing the filling into smaller segments, fragments of the filling might already detach themselves from cavity.
- Filling residues are removed with suitable hand instruments.

5. Reprocessing

Those medical devices that are supplied non-sterile have to be prepared prior to first use. Information provided by the manufacturer regarding reprocessing in compliance with ISO 17664 "Reprocessing of medical devices semi-critical A and B", the present instructions for use and more detailed information, as amended, can be accessed for download under <http://ifr.drendel.com> or requested from the manufacturer Drendel + Zweiling DIAMANT GmbH.

- Initial processing at the site of use/preparation: After use, the product should be placed in a washer/disinfector filled with a suitable detergent/disinfectant solution (e.g. Cidex OPA, co. ASP [validated procedure: undiluted]) to prevent protein fixation. Please observe the manufacturer's information regarding the concentration and temperature as well as minimum and maximum immersion time. After use, reprocessing should take place in a timely manner (≤ 2 hours). The instruments should be inside the washer/disinfector during transport to the reprocessing site. Remove the product from the instrument bath immediately before the start of the cleaning process. Rinse off visible contamination and cleaning/disinfectant agent residues using a nylon brush (do not use a wire brush) under running water (less than 45 °C to prevent protein fixation) for 1 minute.
- Manual cleaning/disinfection in an ultrasonic device: Place the product inside a suitable sieve in an ultrasonic device filled with detergent/disinfectant (e.g. Cidezyme, co. ASP [validated procedure: 0.8 % solution]). Please observe the manufacturer's information regarding the suitability, immersion time, concentration and durability of the solution for ultrasonic cleaning. The instruments must come in contact with each other during the cleaning process in the ultrasonic bath. Thoroughly rinse the instrument with water for 1 minute.
- This is followed by disinfection in the ultrasonic device using a suitable disinfectant (e.g. Cidex OPA, co. ASP [validated procedure: undiluted]).

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Instructions for use

Please read carefully and retain for future reference.

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Make sure to observe the instructions provided by the manufacturer concerning the concentration and immersion time during cleaning and chemical disinfection. The immersion time does not start until the last instrument was placed into the ultrasonic bath. Never exceed the indicated immersion time. Attention: Do not exceed 45 °C (risk of protein coagulation).

- Thoroughly rinse off disinfectant with water and dry carefully.
 - Mechanical cleaning/disinfection in the washer/disinfector: The device must be tested in accordance with ISO 15883. Place the product in the washer/disinfector in a horizontal position, so that the spray jet directly reaches the instruments through the openings. Start the Vario TD program (co. Miele) or the universal program (co. Melag). Observe the manufacturer's information regarding the concentration and quantity of the cleaning agents (e.g. Neodisher MediClean forte, co. Dr. Weigert). Observe the instructions of the manufacturer. Thermal disinfection takes place in the washer/disinfector (at least 5 minutes at 90 °C or A, value ≥ 3000). On completion of the complete cycle, allow instruments to cool down and remove from the washer/disinfector. Dry product carefully.
 - Control/sterilization: Subject the product to a visual control. Do not reuse corroded, damaged or blunt instruments. The sterilization takes place in the autoclave type B (EN 13060/EN 285 or ANSI AAMI ST79/approved by the FDA), tested in accordance with ISO 17665 in suitable sterilization packaging at 134 °C for 5 minutes, and a subsequent drying time of 20 minutes. The information provided by the manufacturer of the device have to be observed. Please also observe the legal provision valid in your country.
- The effectiveness of the described procedure was confirmed by an accredited laboratory. If the user is using a different method, the user has to validate the efficiency of such a method. The operator of a medical device is responsible for ensuring that reprocessing takes place with suitable equipment, suitable materials and by duly qualified staff and that this is carried out in compliance with any legal provisions valid in his or her country.

6. Storage

Store wrapped devices protected from UV radiation and high temperatures in a clean and dry environment. Do not store in the same room as solvents or chemicals.

7. Disposal

The devices have to be disposed of in tight, non-breakable and puncture-proof containers (protection from contamination).

8. Safety and liability

Worn and damaged instruments (bent instruments etc.) have to be rejected and replaced by new ones. The above mentioned recommendations with respect to handling, cooling and contact pressure are to be strictly observed. The instruments must only be used for the intended application. Non-observance of these safety recommendations may lead to damage to the handpiece or injury. Please also observe the General Instructions for Use and Safety Recommendations in the current Drendel+Zweiling catalogue. The user is obliged to check the devices on his own responsibility prior to use to ensure that they are suited for the intended purpose. The correct application of the devices is the responsibility of the operator and user. In case of contributory negligence by the user, Drendel+Zweiling DIAMANT GmbH partially or totally declines liability towards the user for all resulting damages, particularly if these are due to the non-observance of the recommendations for use or other misuse by the user. Store out of children's reach.

9. Contains hazardous substances

Devices marked with a CAS number (CAS: 7440-48-4) contain more than 0.1 % weight percent of cobalt. As a CMR substance of class 1B, cobalt is classified as possibly carcinogenic, mutagenic and/or toxic to reproduction. Evidence has shown that the quantities of cobalt released by medical devices are so low that they do not pose a risk and no precautions have to be taken, provided that the medical device is used correctly and according to its intended purpose.

10. Procedure in case of serious adverse events

In case of serious adverse events, please proceed in line with the legal regulations in your country.