Füllungsbearbeitung, Work on fillings, Preparation des obturations, Trabajo en obturaciones, Lavorazione di otturazioni

Verpackungseinheit, Packing unit, Unité d'emballage, Unidad de embalaje, Unità d'imballo

Enthält gefährliche Substanzen, Contains hazardous substances, Contient des substances dangereuses, Contiene sustancias peligrosas, Contiene sostanze pericolose

7440-48-4: Cobalt, 7440-48-4: Cobalt, 7440-48-4: Cobalt, 7440-48-4: Cobalt, 7440-48-4: Cobalt


Fig. 1


Fig. 2

Kombinationsinstrumente aus Hartmetall

1. Indikation

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Tungsten carbide combination instruments

1. Indication

The instruments are intended for trimming and finishing plastic restoration materials.

2. Contraindication

No product-related contraindications are known.

3. Proper use

- Use instruments after the cavity has been filled with plastic restoration materials.

- Trimming and finishing is carried out step by step.

a) Trim the restoration surface with the selected instrument in a red contra-angle at 20,000 rpm with usual coolant supply (50 ml/min.). The combination finishing instrument is used like a traditional finishing instrument, i.e. cutting to the right, and the contra-angle is guided contrary to the rotation direction of the instrument (Fig. 1: guiding the instrument contrary to its rotary movement). Prepare the filling so that only a small amount of material has to be removed during the subsequent finishing.

b) Use the same instrument for final finishing of the filling in an anticlockwise direction and guide the contra-angle in the rotary movement of the instrument (Fig. 2: guiding the instrument in its rotary movement). Set the motor to anti-clockwise rotation. The recommended speed is 20,000 rpm in a red contra-angle with usual coolant supply (50 ml/min.).

- Apart from this general information, please observe the individual speed indicated on the packaging.

- The maximum permissible speed indicated on the package must not be exceeded in order to avoid risk of injury due to instrument fracture.

4. Reprocessing

Those medical devices that are supplied non-sterile have to be prepared prior to first use. Information provided by the manufacturer regarding reprocessing in compliance with ISO 17664 "Reprocessing of medical devices semi-critical A and B", the present instructions for use and more detailed information, as amended, can be accessed for download under <http://ifr.drendel.com> or requested from the manufacturer Drendel + Zweiling DIAMANT GmbH.

- Initial processing at the site of use/preparation: After use, the product should be placed in a washer/disinfector filled with a suitable detergent/disinfectant solution (e.g. Cidex OPA, co. ASP [validated procedure: undiluted]) to prevent protein fixation. Please observe the manufacturer's information regarding the concentration and temperature as well as minimum and maximum immersion time. After use, reprocessing should take place in a timely manner (≤ 2 hours). The instruments should be inside the washer/disinfector during transport to the reprocessing site. Remove the product from the instrument bath immediately before the start of the cleaning process. Rinse off visible contamination and cleaning/disinfectant agent residues using a nylon brush (do not use a wire brush) under running water (less than 45 °C to prevent protein fixation) for 1 minute.

- Manual cleaning/disinfection in an ultrasonic device: Place the product inside a suitable sieve in an ultrasonic device filled with detergent/disinfectant (e.g. Cidezyme, co. ASP [validated procedure: 0.8 % solution]). Please observe the manufacturer's information regarding the suitability, immersion time, concentration and durability of the solution for ultrasonic cleaning. The instruments must come in contact with each other during the cleaning process in the ultrasonic bath. Thoroughly rinse the instrument with water for 1 minute.

This is followed by disinfection in the ultrasonic device using a suitable disinfectant (e.g. Cidex OPA, co. ASP [validated procedure: undiluted]).

Make sure to observe the instructions provided by the manufacturer concerning the concentration and immersion time during cleaning and chemical disinfection. The immersion time does not start until the last instrument was placed into the ultrasonic bath. Never exceed the

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Instructions for use

Please read carefully and retain for future reference.

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indicated immersion time. Attention: Do not exceed 45 °C (risk of protein coagulation).

Thoroughly rinse off disinfectant with water and dry carefully.

- Mechanical cleaning/disinfection in the washer/disinfector: The device must be tested in accordance with ISO 15883. Place the product in the washer/disinfector in a horizontal position, so that the spray jet directly reaches the instruments through the openings. Start the Vario TD program (co. Miele) or the universal program (co. Melag). Observe the manufacturer's information regarding the concentration and quantity of the cleaning agents (e.g. Neodisher MediClean forte, co. Dr. Weigert). Observe the instructions of the manufacturer. Thermal disinfection takes place in the washer/disinfector (at least 5 minutes at 90 °C or A, value ≥ 3000). On completion of the complete cycle, allow instruments to cool down and remove from the washer/disinfector. Dry product carefully.

- Control/sterilization: Subject the product to a visual control. Do not reuse corroded, damaged or blunt instruments. The sterilization takes place in the autoclave type B (EN 13060/EN 285 or ANSI AAMI ST79/approved by the FDA), tested in accordance with ISO 17665 in suitable sterilization packaging at 134 °C for 5 minutes, and a subsequent drying time of 20 minutes. The information provided by the manufacturer of the device have to be observed. Please also observe the legal provision valid in your country.

The effectiveness of the described procedure was confirmed by an accredited laboratory. If the user is using a different method, the user has to validate the efficiency of such a method. The operator of a medical device is responsible for ensuring that reprocessing takes place with suitable equipment, suitable materials and by duly qualified staff and that this is carried out in compliance with any legal provisions valid in his or her country.

5. Storage

Store wrapped devices protected from UV radiation and high temperatures in a clean and dry environment. Do not store in the same room as solvents or chemicals.

6. Disposal

The devices have to be disposed of in tight, non-breakable and puncture-proof containers (protection from contamination).

7. Safety and liability

Worn and damaged instruments (bent instruments etc.) have to be rejected and replaced by new ones. The above mentioned recommendations with respect to handling, cooling and contact pressure are to be strictly observed. The instruments must only be used for the intended application. Non-observance of these safety recommendations may lead to damage to the handpiece or injury. Please also observe the General Instructions for Use and Safety Recommendations in the current Drendel+Zweiling catalogue. The user is obliged to check the devices on his own responsibility prior to use to ensure that they are suited for the intended purpose. The correct application of the devices is the responsibility of the operator and user. In case of contributory negligence by the user, Drendel+Zweiling DIAMANT GmbH partially or totally declines liability towards the user for all resulting damages, particularly if these are due to the non-observance of the recommendations for use or other misuse by the user. Store out of children's reach. For dental use only.

8. Contains hazardous substances

Devices marked with a CAS number (CAS: 7440-48-4) contain more than 0.1 % weight percent of cobalt. As a CMR substance of class 1B, cobalt is classed as possibly carcinogenic, mutagenic and/or toxic to reproduction. Evidence has shown that the quantities of cobalt released by medical devices are so low that they do not pose a risk and no precautions have to be taken, provided that the medical device is used correctly and according to its intended purpose.

9. Procedure in case of serious adverse events

In case of serious adverse events, please proceed in line with the legal regulations in your country.