ZD Diamond abrasives ZD-Fraises diamantées

ZrO₂ – Diamond abrasives for work on zirconium oxide

The popularity of tooth coloured restorations is constantly on the increase. Zirconium oxide has proven to be the perfect material for this application due to its reliable and durable properties, however, the great hardness of this material makes it very difficult to grind with conventional instruments.

Special ZD diamond abrasives are now available for easy and efficient work on zirconium oxide restorations. Thanks to a special bonding and specially chosen diamond grains, these abrasive instruments facilitate efficient material reduction, while creating a perfect surface finish without dark friction marks.

Special diamond coating:

- Optimal material reduction
- Considerably longer service life compared to conventional diamond instruments

Different shapes are available for:

- Quick trepanation of ZrO₂ restorations with the ZD801
- Minor corrections of crowns using the ZD850
- Removing all-ceramic crowns and bridges with the ZD880CC

Recommended speed: opt. 160.000 rpm

We recommend using these abrasives in the red contra-angle (rather than in the conventional turbine) as the higher torque allows more efficient work on ZrO₂. Make sure to supply sufficient coolant, particularly during trepanation (min. 50 ml/min.).

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DIAMANT GmbH