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ZD-Diamonds
ZD-Instruments
diamantés


Gates Glidden
Gates Glidden

Polishing discs
Disques à polir

Milling technique
Technique de
fraisage

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D+Z offers a complete range of instruments for the dental office and laboratory.

Why look through a lot of catalogues from different suppliers when you can get all you need from one source?

The complete D+Z range assures an optimal basis for safe working results.


Drendel+Zweiling has always been customer oriented and therefore keeps extending its range for the dental office and laboratory constantly.

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ZD Diamond abrasives
ZD Instruments diamantés

**ZrO₂ – Diamond abrasives
 for work on zirconium oxide**

The popularity of tooth coloured restorations is constantly on the increase. Zirconium oxide has proven to be the perfect material for this application due to its reliable and durable properties, however, the great hardness of this material makes it very difficult to grind with conventional instruments.

Thanks to a special bonding and specially chosen diamond grains, these abrasive instruments facilitate efficient material reduction while creating a perfect surface finish without dark friction marks.

For optimal results we recommend to use the instruments in the red contra-angle at an optimal speed of 160,000 rpm.

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Gates Glidden
Gates Glidden

Gates Glidden


Once the root canals have been retrieved, the enlargement of the straight, coronal part of the canal entrance is a crucial step to help all subsequently used instruments penetrate and efficiently prepare the deeper portions of the canal.

Our stainless steel Gates Glidden reamers are applied in the Crown-Down technique and facilitate tapered reaming in the straight, coronal part of the root canal.

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

Composite Polishing *Polissage des composites*

Single-use polishers for composite restorations

The polishing discs and polishing strips facilitate optimal polishing of composite in three steps.


The transparent discs are provided with an aluminium oxide coating on both sides. The discs are mounted onto the Pop-on mandrel to assure precise and secure locking. The double-sided coating of the discs allows pushing and pulling movements without time-consuming change-over. The mandrel is completely covered by the discs, which minimizes the risk of damaging the restoration.

The strips are based on a clever 3 in 1 principle. All three grit sizes (from fine to coarse) are arranged in succession on each strip. Uncoated surfaces for insertion and handy grips at the ends facilitate the use of these polishing strips.

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Milling technique Technique de fraisage

Milling technique

A new range of cutters and polishers is now available from D+Z to facilitate efficient work on parallel and tapered surfaces (such as telescopic crowns, tapered crowns, bars for implant suprastructures and abutments) using the milling device.

Wax cutters CM364RA and CM356RA

The new wax cutters CM364RA and CM356RA achieve very fine surfaces in no time at all. Thanks to their special spiral-shaped blade geometry, subsequent scraping can be omitted. Important: Synchronous milling is recommended (in clockwise direction and with a clockwise drive) at a reduced speed of 3,000 rpm.

Coarse cutter CM364RCX and CM356RCX for non-precious metal alloys and titanium

The cutters CM364RCX and CM356RCX are designed for initial trimming of hard alloys. The cutters feature a coarse staggered toothing for optimal material reduction. Thanks to their prolonged service life they are economic in use. The cutters are operated in a contra-rotational direction using a clockwise drive. You can reduce premature wear of your cutter by using milling oil and working at a reduced speed of 6,000 rpm.

Fine cutters CM364RS and CM356RS for non-precious metal alloys and titanium

Used subsequently to coarse cutters, the fine cutters CM364RS and CM365RS achieve shiny surfaces thanks to their efficient toothing. The cutters are used in a contra-rotational direction at a reduced speed of 3,000 rpm and have to be lubricated with oil. If required, rework using oil and wax.

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





Coarse cutters CM364RMX and CM356RMX for gold alloys

D+Z recommends the cutters CM364RMX and CM356RMX with medium staggered toothing for the pre-milling of soft alloys or gold alloys. The cutters are also used in a contra-rotational direction at an optimal speed of 10,000 rpm and have to be lubricated with oil.

Fine cutters CM364RF and CM356RF

D+Z recommends the cutters CM364RF and CM356RF for fine trimming or initial polishing. The toothing of these cutters is provided with a special chamfer and is ideally suited for work on gold alloys. The cutters are used in contra-rotational direction at a speed of 3,000 rpm and have to be lubricated with oil and wax, if required.

Polishers for milling technique P9440C/M/F

The polishers P9440C, P9440M and P9440F are designed for polishing telescopic crowns. Prior to polishing, the surfaces of crowns have to be finely milled to ensure that they are smooth and scratch-free. Polishing is carried out in three steps (C = coarse (dark brown), M = medium (reddish brown), F = fine (green)) from initial polishing to high-shine polishing. The polishers are used dry at a maximum speed of 6,000 rpm.

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Separating disc 9507 Disque à séparer 9507

Fibre-reinforced separating disc 9507

This separating disc is predestined for efficient and quick removal of sprues on non-precious metal alloys. The disc is reinforced with glass fibre which makes it particularly resistant to fracture.

The disc in size 250 is used at an optimal speed of 20,000 rpm, whereas size 400 is preferably used at a reduced speed of 10,000 rpm.

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DIAMANT GmbH