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EN

Instructions for use and safety recommendations for the application of rotary instruments in the medical field

D+Z instruments are made of high-grade materials. They are designed, constructed and produced with utmost care. These quality instruments assure best working results and long service life provided they are used in the proper manner. Therefore, the following instructions for use and safety recommendations have to be observed. Improper use can lead to damage to the tissue, premature wear, destruction of the instruments and injury to the operator, patient or other persons.

1. Proper use

- Observe this information and the instructions for use for surgical motor systems when using rotary instruments.
- Prior to use, make sure that the instruments are in perfect technical condition and that sterility is assured. Damaged instruments should not be reused.
- Insert rotary instruments in the hand piece until the stop position is reached.
- Instruments with clamping marks (line or ring) have to be lowered into the hand-piece until the marks are no longer visible.
- Prior to operation, make sure that the instrument is firmly locked.
- Do not use the instruments for work on metal (e.g. implants) and do not allow them to come in contact with metal during use as this might cause nicks in the blades.
- Craniotomes are to be used with corresponding dura guard only.
- To ensure traceability of the MD by means of a UDI, please make sure to keep the packaging of articles without direct marking.

2. Indication

D+Z instruments have been designed for working on bone, skin, soft tissue and cements. Instruments of different materials adapted to each application are available.
- Maxillo-facial surgery:
instruments made of stainless steel, tungsten carbide instruments, diamond instruments
- ENT surgery:
instruments made of stainless steel, tungsten carbide instruments, diamond instruments
- Neurosurgery:
stainless steel, carbide and diamond instruments
- Orthopedics:
stainless steel, carbide and diamond instruments
The surgeon is responsible for the selection of the appropriate instrument.

3. Contraindication

The instruments are not suited for use in the CNS (brain and spinal cord).

4. Contact pressure

- Excessive contact pressure has to be avoided.
- It can lead to damage of the instruments' working part (steel and tungsten carbide instruments) as well as chipping of the . At the same time extreme heat is generated which increases the risk of thermal necrosis.
- Excessive contact pressure can lead to the abrasive grit of diamond instruments breaking off. Increased heat development will be the result.
- At worst, breakage of the instrument can occur if excessive contact pressure is applied. Jamming or leverage of the instrument must be avoided during operation.
- Optimal material reduction is assured at a contact pressure of approx. 2N. This can be compared with the pressure applied when writing with a fountain pen.

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Instructions for use

Please read carefully and retain for future reference.

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5. Cooling

- When using rotary instruments sufficient coolinghas to be supplied, otherwise the instrument blades will be clogged with chips which leads to increased heat generation and at worst to irreversible damage to the bone (thermal necrosis). Moreover, the service life of the instruments will be reduced.

6. Elimination of blunt and damaged instruments

- Examine instruments for bluntness and damage with magnifying glasses. Look for the following indications:
- blunt (rounded) blades and chipping
 - blank spots on diamond instruments
 - damage to the shank
 - bent or "out of true" instruments
- Do not reuse blunt or damaged instruments under any circumstances. This might lead to fracture of the instrument and risk of injury to the patient and the operator.

7. Cleaning, disinfection, sterilization and storage

- A validated reprocessing method is available on our web site www.drendel.com!
 - Prior to first use and immediately after each use, rotary instruments have to be disinfected, cleaned and sterilized. The instruments should be stored in their original packages at room temperature, protected against dust and humidity.
 - The instruments should be placed in hygienic bur blocks, dishes or other suitable containers. The same applies to sterilized instruments and instruments in sterile packages. The instruments are to be stored in an environment free of dust and humidity, to exclude recontamination.
 - After surgical operation, immediately remove blood, secretion, tissue and bone residues. Dried residues can cause corrosion.
- Make sure that all residues are removed; this applies especially to hollow burs.
- D+Z instruments can be cleaned in the ultrasonic bath. However, care should be taken that the instruments do not come in contact with each other; otherwise the blades can be damaged.
 - It is recommended to use anticorrosive disinfection and cleaning agents.
 - Contact with H₂O₂ (hydrogen peroxide), chlorous disinfection and cleaning agents as well as oxalic acids are to be avoided because these agents may lead to pitting and corrosion.
 - For immersion time and suitability of disinfection and cleaning substances for specific instrument types see instructions for use of the manufacturer of these agents.
 - Disinfection and cleaning agents must be rinsed off thoroughly under running water immediately after use; otherwise corrosion can occur due to chemical reactions.
 - Disinfection of rotary instruments can be carried out in the gentle cycle (e.g. Vario TD) of the thermodesinfector. Use cleaning and disinfection agents with approved material compatibility. We recommend deconex 28 AlkaOne (co. Borer Chemie).
 - Sterilization should take place in the steam sterilizer.
 - A temperature above 180°C has to be avoided. Exceeding this temperature would result in loss of hardness of the working part and thus reduction of service life.
- Make sure that each selected disinfection and sterilization method is suitable for the specific instruments. The instruments may only be returned in clean, disinfected and sterilized condition.

8. Safety and liability

- The user is responsible to check the product prior to use whether it is suitable for the intended purpose.
- The user is responsible for the application of the instruments.
- In case of contributory negligence by the user, D+Z partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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© 1218.304327 – Rev 12/2018